

RT² Profiler PCR Array (Rotor-Gene[®] Format)

Human Extracellular Matrix & Adhesion Molecules

Cat. no. 330231 PAHS-013ZR

For pathway expression analysis

Format	For use with the following real-time cyclers
RT ² Profiler PCR Array, Format R	Rotor-Gene Q, other Rotor-Gene cyclers

Description

The Human Extracellular Matrix & Adhesion Molecules RT² Profiler PCR Array profiles the expression of 84 genes important for cell-cell and cell-matrix interactions. Cells attach to proteoglycans and glycoproteins (such as fibronectin, laminin, and collagen) in the extracellular matrix (ECM) substratum via adhesion molecules on their cell surface to define tissue shape, structure, and function. Making and breaking cellular contacts with other cells and the ECM play critical roles in normal processes such as cell growth, division, differentiation, and migration. Disease states (metastasis, rheumatoid arthritis, and cardiovascular and central nervous system disorders) and pathophysiological processes (wound healing and inflammation) also involve ECM remodeling and changes in the classes of cell adhesion molecule expressed on cell surfaces. This array contains ECM proteins including basement membrane constituents, collagens, and other genes defining ECM structure. Matrix and other metalloproteinases that remodel the ECM as well as their inhibitors are also included. This array also represents integrins, selectins, cell-adhesion molecule family members (ICAM, ECAM, NCAM, PECAM, and VCAM), and other genes important to cell adhesion and cytoskeleton bridging such as the catenins. Using real-time PCR, you can easily and reliably analyze expression of a focused panel of genes related to cell adhesion and the ECM with this array.

For further details, consult the *RT² Profiler PCR Array Handbook*.


Shipping and storage

RT² Profiler PCR Arrays in the Rotor-Gene format are shipped at ambient temperature, on dry ice, or blue ice packs depending on destination and accompanying products.

For long term storage, keep plates at -20°C .

Note: Ensure that you have the correct RT² Profiler PCR Array format for your real-time cyclers (see table above).

Note: Open the package and store the products appropriately immediately on receipt.



Array layout

The 96 real-time assays in the Rotor-Gene format are located in wells 1–96 of the Rotor-Disc™ (plate A1–A12=Rotor-Disc 1–12, plate B1–B12=Rotor-Disc 13–24, etc.). To maintain data analysis compatibility, wells 97–100 do not contain real-time assays but will contain master mix to account for weight balance.

Gene table: RT² Profiler PCR Array

Position	UniGene	GenBank	Symbol	Description
A01	Hs.643357	NM_006988	ADAMTS1	ADAM metalloproteinase with thrombospondin type 1 motif, 1
A02	Hs.131433	NM_139025	ADAMTS13	ADAM metalloproteinase with thrombospondin type 1 motif, 13
A03	Hs.271605	NM_007037	ADAMTS8	ADAM metalloproteinase with thrombospondin type 1 motif, 8
A04	Hs.502328	NM_000610	CD44	CD44 molecule (Indian blood group)
A05	Hs.461086	NM_004360	CDH1	Cadherin 1, type 1, E-cadherin (epithelial)
A06	Hs.476092	NM_003278	CLEC3B	C-type lectin domain family 3, member B
A07	Hs.143434	NM_001843	CNTN1	Contactin 1
A08	Hs.523446	NM_080629	COL11A1	Collagen, type XI, alpha 1
A09	Hs.101302	NM_004370	COL12A1	Collagen, type XII, alpha 1
A10	Hs.409662	NM_021110	COL14A1	Collagen, type XIV, alpha 1
A11	Hs.409034	NM_001855	COL15A1	Collagen, type XV, alpha 1
A12	Hs.368921	NM_001856	COL16A1	Collagen, type XVI, alpha 1
B01	Hs.172928	NM_000088	COL1A1	Collagen, type I, alpha 1
B02	Hs.508716	NM_001846	COL4A2	Collagen, type IV, alpha 2
B03	Hs.210283	NM_000093	COL5A1	Collagen, type V, alpha 1
B04	Hs.474053	NM_001848	COL6A1	Collagen, type VI, alpha 1
B05	Hs.420269	NM_001849	COL6A2	Collagen, type VI, alpha 2
B06	Hs.476218	NM_000094	COL7A1	Collagen, type VII, alpha 1
B07	Hs.654548	NM_001850	COL8A1	Collagen, type VIII, alpha 1
B08	Hs.591346	NM_001901	CTGF	Connective tissue growth factor
B09	Hs.534797	NM_001903	CTNNA1	Catenin (cadherin-associated protein), alpha 1, 102kDa
B10	Hs.476018	NM_001904	CTNNB1	Catenin (cadherin-associated protein), beta 1, 88kDa
B11	Hs.166011	NM_001331	CTNND1	Catenin (cadherin-associated protein), delta 1
B12	Hs.314543	NM_001332	CTNND2	Catenin (cadherin-associated protein), delta 2 (neural plakophilin-related arm-repeat protein)
C01	Hs.81071	NM_004425	ECM1	Extracellular matrix protein 1
C02	Hs.203717	NM_002026	FN1	Fibronectin 1
C03	Hs.57697	NM_001523	HAS1	Hyaluronan synthase 1
C04	Hs.643447	NM_000201	ICAM1	Intercellular adhesion molecule 1
C05	Hs.644352	NM_181501	ITGA1	Integrin, alpha 1
C06	Hs.482077	NM_002203	ITGA2	Integrin, alpha 2 (CD49B, alpha 2 subunit of VLA-2 receptor)
C07	Hs.265829	NM_002204	ITGA3	Integrin, alpha 3 (antigen CD49C, alpha 3 subunit of VLA-3 receptor)
C08	Hs.694732	NM_000885	ITGA4	Integrin, alpha 4 (antigen CD49D, alpha 4 subunit of VLA-4 receptor)
C09	Hs.505654	NM_002205	ITGA5	Integrin, alpha 5 (fibronectin receptor, alpha polypeptide)
C10	Hs.133397	NM_000210	ITGA6	Integrin, alpha 6
C11	Hs.524484	NM_002206	ITGA7	Integrin, alpha 7
C12	Hs.171311	NM_003638	ITGA8	Integrin, alpha 8
D01	Hs.174103	NM_002209	ITGAL	Integrin, alpha L (antigen CD11A (p180), lymphocyte function-associated antigen 1; alpha polypeptide)
D02	Hs.172631	NM_000632	ITGAM	Integrin, alpha M (complement component 3 receptor 3 subunit)
D03	Hs.436873	NM_002210	ITGAV	Integrin, alpha V (vitronectin receptor, alpha polypeptide, antigen CD51)
D04	Hs.643813	NM_002211	ITGB1	Integrin, beta 1 (fibronectin receptor, beta polypeptide, antigen CD29 includes MDF2, MSK12)
D05	Hs.375957	NM_000211	ITGB2	Integrin, beta 2 (complement component 3 receptor 3 and 4 subunit)
D06	Hs.218040	NM_000212	ITGB3	Integrin, beta 3 (platelet glycoprotein IIIa, antigen CD61)
D07	Hs.632226	NM_000213	ITGB4	Integrin, beta 4
D08	Hs.536663	NM_002213	ITGB5	Integrin, beta 5
D09	Hs.521869	NM_000216	KAL1	Kallmann syndrome 1 sequence
D10	Hs.270364	NM_005559	LAMA1	Laminin, alpha 1
D11	Hs.200841	NM_000426	LAMA2	Laminin, alpha 2
D12	Hs.436367	NM_000227	LAMA3	Laminin, alpha 3
E01	Hs.650585	NM_002291	LAMB1	Laminin, beta 1
E02	Hs.497636	NM_000228	LAMB3	Laminin, beta 3
E03	Hs.609663	NM_002293	LAMC1	Laminin, gamma 1 (formerly LAMB2)
E04	Hs.83169	NM_002421	MMP1	Matrix metalloproteinase 1 (interstitial collagenase)
E05	Hs.2258	NM_002425	MMP10	Matrix metalloproteinase 10 (stromelysin 2)
E06	Hs.143751	NM_005940	MMP11	Matrix metalloproteinase 11 (stromelysin 3)

Position	UniGene	GenBank	Symbol	Description
E07	Hs.1695	NM_002426	MMP12	Matrix metalloproteinase 12 (macrophage elastase)
E08	Hs.2936	NM_002427	MMP13	Matrix metalloproteinase 13 (collagenase 3)
E09	Hs.2399	NM_004995	MMP14	Matrix metalloproteinase 14 (membrane-inserted)
E10	Hs.80343	NM_002428	MMP15	Matrix metalloproteinase 15 (membrane-inserted)
E11	Hs.546267	NM_005941	MMP16	Matrix metalloproteinase 16 (membrane-inserted)
E12	Hs.513617	NM_004530	MMP2	Matrix metalloproteinase 2 (gelatinase A, 72kDa gelatinase, 72kDa type IV collagenase)
F01	Hs.375129	NM_002422	MMP3	Matrix metalloproteinase 3 (stromelysin 1, progelatinase)
F02	Hs.2256	NM_002423	MMP7	Matrix metalloproteinase 7 (matrilysin, uterine)
F03	Hs.161839	NM_002424	MMP8	Matrix metalloproteinase 8 (neutrophil collagenase)
F04	Hs.297413	NM_004994	MMP9	Matrix metalloproteinase 9 (gelatinase B, 92kDa gelatinase, 92kDa type IV collagenase)
F05	Hs.503878	NM_000615	NCAM1	Neural cell adhesion molecule 1
F06	Hs.514412	NM_000442	PECAM1	Platelet/endothelial cell adhesion molecule
F07	Hs.89546	NM_000450	SELE	Selectin E
F08	Hs.728756	NM_000655	SELL	Selectin L
F09	Hs.73800	NM_003005	SELP	Selectin P (granule membrane protein 140kDa, antigen CD62)
F10	Hs.371199	NM_003919	SGCE	Sarcoglycan, epsilon
F11	Hs.111779	NM_003118	SPARC	Secreted protein, acidic, cysteine-rich (osteonectin)
F12	Hs.185597	NM_003119	SPG7	Spastic paraplegia 7 (pure and complicated autosomal recessive)
G01	Hs.313	NM_000582	SPP1	Secreted phosphoprotein 1
G02	Hs.369397	NM_000358	TGFB1	Transforming growth factor, beta-induced, 68kDa
G03	Hs.164226	NM_003246	THBS1	Thrombospondin 1
G04	Hs.371147	NM_003247	THBS2	Thrombospondin 2
G05	Hs.169875	NM_007112	THBS3	Thrombospondin 3
G06	Hs.522632	NM_003254	TIMP1	TIMP metalloproteinase inhibitor 1
G07	Hs.633514	NM_003255	TIMP2	TIMP metalloproteinase inhibitor 2
G08	Hs.644633	NM_000362	TIMP3	TIMP metalloproteinase inhibitor 3
G09	Hs.143250	NM_002160	TNC	Tenascin C
G10	Hs.109225	NM_001078	VCAM1	Vascular cell adhesion molecule 1
G11	Hs.643801	NM_004385	VCAN	Versican
G12	Hs.2257	NM_000638	VTN	Vitronectin
H01	Hs.520640	NM_001101	ACTB	Actin, beta
H02	Hs.534255	NM_004048	B2M	Beta-2-microglobulin
H03	Hs.592355	NM_002046	GAPDH	Glyceraldehyde-3-phosphate dehydrogenase
H04	Hs.412707	NM_000194	HPRT1	Hypoxanthine phosphoribosyltransferase 1
H05	Hs.546285	NM_001002	RPLP0	Ribosomal protein, large, P0
H06	N/A	SA_00105	HGDC	Human Genomic DNA Contamination
H07	N/A	SA_00104	RTC	Reverse Transcription Control
H08	N/A	SA_00104	RTC	Reverse Transcription Control
H09	N/A	SA_00104	RTC	Reverse Transcription Control
H10	N/A	SA_00103	PPC	Positive PCR Control
H11	N/A	SA_00103	PPC	Positive PCR Control
H12	N/A	SA_00103	PPC	Positive PCR Control

Related products

For optimal performance, RT² Profiler PCR Arrays should be used together with the RT² First Strand Kit for cDNA synthesis and RT² SYBR[®] Green qPCR Mastermixes for PCR.

Product	Contents	Cat. no.
RT ² First Strand Kit (12)	Enzymes and reagents for cDNA synthesis	330401
RT ² SYBR Green ROX™ FAST Mastermix (2)*	For 2 x 96 assays in 96-well plates; suitable for use with the Rotor-Gene Q and other Rotor-Gene cyclers	330620

* Larger kit sizes available; please inquire.

RT² Profiler PCR Array products are intended for molecular biology applications. These products are not intended for the diagnosis, prevention, or treatment of a disease.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN[®], Rotor-Gene[®], Rotor-Disc[™] (QIAGEN Group); ROX[™] (Applied Biosystems Corporation or its subsidiaries); SYBR[®] (Molecular Probes, Inc.).

1067688 03/2011 © 2011 QIAGEN, all rights reserved.

www.qiagen.com

Australia ■ 1-800-243-800

Austria ■ 0800/281010

Belgium ■ 0800-79612

Brazil ■ 0800-557779

Canada ■ 800-572-9613

China ■ 8621-3865-3865

Denmark ■ 80-885945

Finland ■ 0800-914416

France ■ 01-60-920-930

Germany ■ 02103-29-12000

Hong Kong ■ 800 933 965

Ireland ■ 1800 555 049

Italy ■ 800-787980

Japan ■ 03-6890-7300

Korea (South) ■ 080-000-7145

Luxembourg ■ 8002 2076

Mexico ■ 01-800-7742-436

The Netherlands ■ 0800 0229592

Norway ■ 800-18859

Singapore ■ 1800-742-4368

Spain ■ 91-630-7050

Sweden ■ 020-790282

Switzerland ■ 055-254-22-11

UK ■ 01293-422-911

USA ■ 800-426-8157



Sample & Assay Technologies