

QuantiNova® LNA® Probe PCR Focus Panels (Rotor-Gene® Format)

Human Th1 & Th2 Responses

Cat. no. 249955 UPHS-034ZR

For study focus gene expression analysis

Shipping and storage

QuantiNova LNA Probe PCR Focus Panels are shipped at room temperature. Immediately upon receipt, they should be stored protected from light at 2–8°C for short term storage or at –30°C to –15°C for long time storage. Under these conditions, all components are stable for at least 12 months.

Note: Open the package and store the products appropriately immediately upon receipt.

For optimal performance, QuantiNova LNA Probe PCR Focus Panels should be used together with the QuantiNova Reverse Transcription Kit for cDNA synthesis and the QuantiNova Probe PCR Kit (Mastermix) for PCR.

Panel layout (Rotor-Gene): QuantiNova LNA Probe PCR Focus Panel

The 96 real-time assays in the Rotor-Gene format are located in wells 1–96 of the Rotor-Disc® (plate A1–A12=Rotor-Disc 1–12, plate B1–B12=Rotor-Disc 13–24, etc.). To maintain data analysis compatibility, wells 97–100 do not contain real-time assays but will contain master mix to account for weight balance. Refer to the QuantiNova LNA Probe PCR Handbook at www.qiagen.com for further details.

	1	2	3	4	5	6	7	8	9	10	11	12
A	BCL6	CCL11	CCL5	CCL7	CCR2	CCR3	CCR4	CCR5	CD27	CD28	CD4	CD40LG
B	CD80	CD86	CEBPB	CREBBP	CSF2	CTLA4	CXCR3	EBI3	FASLG	GATA3	GFI1	PTGDR2
C	HAVCR2	ICOS	IFNG	IL10	IL12B	IL12RB2	IL13	IL13RA1	IL15	IL18	IL18R1	IL1R1
D	IL1RL1	IL2	IL24	IL25	IL27	IL27RA	IL2RA	IL3	IL4	IL4R	IL5	IL6
E	IL6R	IL7	IL7R	IL9	IRF1	IRF4	JAK1	JAK2	LAG3	LAT	LTA	MAF
F	MAPK8	NFATC1	NFATC2	PCGF2	PTPRC	SFTPD	SLC11A1	SOCS1	SOCS5	SPP1	STAT1	STAT4
G	STAT6	TBX21	TGFB3	TLR4	TLR6	TNF	TNFRSF8	TNFRSF9	TNFSF4	TYK2	VEGFA	YY1
H	ACTB	B2M	GAPDH	HPRT1	RPLP0	HGDC	QIC	QIC	QIC	PPC	PPC	PPC

Gene table: QuantiNova LNA Probe PCR Focus Panel

Position	Assay	Name	Symbol	Ensembl ID	Description
A01	UPFH1132778	ENST00000450123.6	BCL6	ENSG00000113916	BCL6, transcription repressor Source HGNC Symbol Acc HGNC 1001
A02	UPFH0201571	ENST00000305869.3	CCL11	ENSG00000172156	C-C motif chemokine ligand 11 Source HGNC Symbol Acc HGNC 10610
A03	UPFH1132786	ENST00000603197.6	CCL5	ENSG00000271503	C-C motif chemokine ligand 5 Source HGNC Symbol Acc HGNC 10632
A04	UPFH0186519	ENST00000378569.2	CCL7	ENSG00000108688	C-C motif chemokine ligand 7 Source HGNC Symbol Acc HGNC 10634
A05	UPFH0175349	ENST00000445132.2	CCR2	ENSG00000121807	C-C motif chemokine receptor 2 Source HGNC Symbol Acc HGNC 1603
A06	UPFH1132788	ENST00000395940.3	CCR3	ENSG00000183625	C-C motif chemokine receptor 3 Source HGNC Symbol Acc HGNC 1604
A07	UPFH0179708	ENST00000330953.5	CCR4	ENSG00000183813	C-C motif chemokine receptor 4 Source HGNC Symbol Acc HGNC 1605
A08	UPFH1132860	ENST00000292303.4	CCR5	ENSG00000160791	C-C motif chemokine receptor 5 (gene/pseudogene) Source HGNC Symbol Acc HGNC 1606
A09	UPFH0539172	ENST00000266557.3	CD27	ENSG00000139193	CD27 molecule Source HGNC Symbol Acc HGNC 11922
A10	UPFH0310921	ENST00000458610.6	CD28	ENSG00000178562	CD28 molecule Source HGNC Symbol Acc HGNC 1653
A11	UPFH1132302	ENST00000541982.5	CD4	ENSG00000101610	CD4 molecule Source HGNC Symbol Acc HGNC 1678
A12	UPFH0592498	ENST00000370629.6	CD40LG	ENSG00000102245	CD40 ligand Source HGNC Symbol Acc HGNC 11935
B01	UPFH1132790	ENST00000264246.8	CD80	ENSG00000121594	CD80 molecule Source HGNC Symbol Acc HGNC 1700
B02	UPFH0045195	ENST00000393627.6	CD86	ENSG00000114013	CD86 molecule Source HGNC Symbol Acc HGNC 1705
B03	UPFH0202295	ENST00000303004.4	CEBPB	ENSG00000172216	CCAAT enhancer binding protein beta Source HGNC Symbol Acc HGNC 1834
B04	UPFH0338543	ENST00000573517.6	CREBBP	ENSG00000005339	CREB binding protein Source HGNC Symbol Acc HGNC 2348
B05	UPFH1132793	ENST00000296871.4	CSF2	ENSG00000164400	colony stimulating factor 2 Source HGNC Symbol Acc HGNC 2434
B06	UPFH0603710	ENST00000427473.3	CTLA4	ENSG00000163599	cytotoxic T-lymphocyte associated protein 4 Source HGNC Symbol Acc HGNC 2505
B07	UPFH1132799	ENST00000373693.4	CXCR3	ENSG00000186810	C-X-C motif chemokine receptor 3 Source HGNC Symbol Acc HGNC 4540
B08	UPFH0245178	ENST00000599339.1	EBI3	ENSG00000105246	Epstein-Barr virus induced 3 Source HGNC Symbol Acc HGNC 3129
B09	UPFH1132396	ENST00000367721.3	FASLG	ENSG00000117560	Fas ligand Source HGNC Symbol Acc HGNC 11936
B10	UPFH1132416	ENST00000645492.1	GATA3	ENSG00000107485	GATA binding protein 3 Source HGNC Symbol Acc HGNC 4172
B11	UPFH0377131	ENST00000294702.6	GFI1	ENSG00000162676	growth factor independent 1 transcriptional repressor Source HGNC Symbol Acc HGNC 4237
B12	UPFH0150140	ENST00000332539.5	PTGDR2	ENSG00000183134	prostaglandin D2 receptor 2 Source HGNC Symbol Acc HGNC 4502
C01	UPFH0231044	ENST00000522593.5	HAVCR2	ENSG00000135077	hepatitis A virus cellular receptor 2 Source HGNC Symbol Acc HGNC 18437
C02	UPFH0092255	ENST00000316386.10	ICOS	ENSG00000163600	inducible T cell costimulator Source HGNC Symbol Acc HGNC 5351
C03	UPFH1132473	ENST00000229135.4	IFNG	ENSG00000111537	interferon gamma Source HGNC Symbol Acc HGNC 5438
C04	UPFH0028177	ENST00000423557.1	IL10	ENSG00000136634	interleukin 10 Source HGNC Symbol Acc HGNC 5962
C05	UPFH0131869	ENST00000231228.2	IL12B	ENSG00000113302	interleukin 12B Source HGNC Symbol Acc HGNC 5970
C06	UPFH0010120	ENST00000541374.5	IL12RB2	ENSG00000081985	interleukin 12 receptor subunit beta 2 Source HGNC Symbol Acc HGNC 5972
C07	UPFH1132807	ENST00000617259.2	IL13	ENSG00000169194	interleukin 13 Source HGNC Symbol Acc HGNC 5973
C08	UPFH0240636	ENST00000371642.1	IL13RA1	ENSG00000131724	interleukin 13 receptor subunit alpha 1 Source HGNC Symbol Acc HGNC 5974
C09	UPFH1132873	ENST00000296545.11	IL15	ENSG00000164136	interleukin 15 Source HGNC Symbol Acc HGNC 5977
C10	UPFH1132481	ENST00000528832.1	IL18	ENSG00000150782	interleukin 18 Source HGNC Symbol Acc HGNC 5986
		ENST00000233		ENSG000000	

Position	Assay	Name	Symbol	Ensembl ID	Description
C11	UPFH0567076	957.5	IL18R1	115604	interleukin 18 receptor 1 Source HGNC Symbol Acc HGNC 5988
C12	UPFH1132482	ENST00000442590.5	IL1R1	ENSG00000115594	interleukin 1 receptor type 1 Source HGNC Symbol Acc HGNC 5993
D01	UPFH0302915	ENST00000404917.6	IL1RL1	ENSG00000115602	interleukin 1 receptor like 1 Source HGNC Symbol Acc HGNC 5998
D02	UPFH0116492	ENST00000226730.4	IL2	ENSG00000109471	interleukin 2 Source HGNC Symbol Acc HGNC 6001
D03	UPFH0437635	ENST00000294984.6	IL24	ENSG00000162892	interleukin 24 Source HGNC Symbol Acc HGNC 11346
D04	UPFH0387713	ENST00000329715.2	IL25	ENSG00000166090	interleukin 25 Source HGNC Symbol Acc HGNC 13765
D05	UPFH0006752	ENST00000356897.1	IL27	ENSG00000197272	interleukin 27 Source HGNC Symbol Acc HGNC 19157
D06	UPFH0315019	ENST00000263379.4	IL27RA	ENSG00000104998	interleukin 27 receptor subunit alpha Source HGNC Symbol Acc HGNC 17290
D07	UPFH0323649	ENST00000447847.1	IL2RA	ENSG00000134460	interleukin 2 receptor subunit alpha Source HGNC Symbol Acc HGNC 6008
D08	UPFH0282899	ENST00000296870.2	IL3	ENSG00000164399	interleukin 3 Source HGNC Symbol Acc HGNC 6011
D09	UPFH0226437	ENST00000231449.7	IL4	ENSG00000113520	interleukin 4 Source HGNC Symbol Acc HGNC 6014
D10	UPFH0363455	ENST00000563886.1	IL4R	ENSG00000077238	interleukin 4 receptor Source HGNC Symbol Acc HGNC 6015
D11	UPFH1132811	ENST00000231454.6	IL5	ENSG00000113525	interleukin 5 Source HGNC Symbol Acc HGNC 6016
D12	UPFH1172910	ENST00000258743.10	IL6	ENSG00000136244	interleukin 6 Source HGNC Symbol Acc HGNC 6018
E01	UPFH1132484	ENST00000622330.4	IL6R	ENSG00000160712	interleukin 6 receptor Source HGNC Symbol Acc HGNC 6019
E02	UPFH1132812	ENST00000263851.9	IL7	ENSG00000104432	interleukin 7 Source HGNC Symbol Acc HGNC 6023
E03	UPFH0314857	ENST00000514217.5	IL7R	ENSG00000168685	interleukin 7 receptor Source HGNC Symbol Acc HGNC 6024
E04	UPFH0611711	ENST00000274520.1	IL9	ENSG00000145839	interleukin 9 Source HGNC Symbol Acc HGNC 6029
E05	UPFH1132490	ENST00000476613.1	IRF1	ENSG00000125347	interferon regulatory factor 1 Source HGNC Symbol Acc HGNC 6116
E06	UPFH0051290	ENST00000380956.9	IRF4	ENSG00000137265	interferon regulatory factor 4 Source HGNC Symbol Acc HGNC 6119
E07	UPFH1132963	ENST00000342505.5	JAK1	ENSG00000162434	Janus kinase 1 Source HGNC Symbol Acc HGNC 6190
E08	UPFH1132818	ENST00000381652.3	JAK2	ENSG00000096968	Janus kinase 2 Source HGNC Symbol Acc HGNC 6192
E09	UPFH0614886	ENST00000203629.3	LAG3	ENSG00000089692	lymphocyte activating 3 Source HGNC Symbol Acc HGNC 6476
E10	UPFH0374827	ENST00000566270.5	LAT	ENSG00000213658	linker for activation of T cells Source HGNC Symbol Acc HGNC 18874
E11	UPFH1132824	ENST00000454783.5	LTA	ENSG00000226979	lymphotxin alpha Source HGNC Symbol Acc HGNC 6709
E12	UPFH0572800	ENST00000326043.5	MAF	ENSG00000178573	MAF bZIP transcription factor Source HGNC Symbol Acc HGNC 6776
F01	UPFH1132535	ENST00000374179.8	MAPK8	ENSG00000107643	mitogen-activated protein kinase 8 Source HGNC Symbol Acc HGNC 6881
F02	UPFH0595445	ENST00000591814.5	NFATC1	ENSG00000131196	nuclear factor of activated T cells 1 Source HGNC Symbol Acc HGNC 7775
F03	UPFH0101128	ENST00000610033.5	NFATC2	ENSG00000101096	nuclear factor of activated T cells 2 Source HGNC Symbol Acc HGNC 7776
F04	UPFH0532144	ENST00000618506.1	PCGF2	ENSG00000277258	polycomb group ring finger 2 Source HGNC Symbol Acc HGNC 12929
F05	UPFH0448301	ENST00000367367.8	PTPRC	ENSG00000081237	protein tyrosine phosphatase, receptor type C Source HGNC Symbol Acc HGNC 9666
F06	UPFH1132678	ENST00000372292.8	SFTPD	ENSG00000133661	surfactant protein D Source HGNC Symbol Acc HGNC 10803
F07	UPFH1132682	ENST00000233202.11	SLC11A1	ENSG00000018280	solute carrier family 11 member 1 Source HGNC Symbol Acc HGNC 10907
F08	UPFH1132887	ENST00000644787.1	SOCS1	ENSG00000185338	suppressor of cytokine signaling 1 Source HGNC Symbol Acc HGNC 19383
F09	UPFH0133368	ENST00000306503.5	SOCS5	ENSG00000171150	suppressor of cytokine signaling 5 Source HGNC Symbol Acc HGNC 16852
F10	UPFH0044238	ENST00000237623.11	SPP1	ENSG00000118785	secreted phosphoprotein 1 Source HGNC Symbol Acc HGNC 11255

Position	Assay	Name	Symbol	Ensembl ID	Description
F11	UPFH1132696	ENST00000392323.6	STAT1	ENSG00000115415	signal transducer and activator of transcription 1 Source HGNC Symbol Acc HGNC 11362
F12	UPFH1132697	ENST00000392320.7	STAT4	ENSG00000138378	signal transducer and activator of transcription 4 Source HGNC Symbol Acc HGNC 11365
G01	UPFH1132845	ENST00000553532.2	STAT6	ENSG00000166888	signal transducer and activator of transcription 6 Source HGNC Symbol Acc HGNC 11368
G02	UPFH0331267	ENST00000177694.2	TBX21	ENSG00000073861	T-box 21 Source HGNC Symbol Acc HGNC 11599
G03	UPFH0000256	ENST00000238682.7	TGFB3	ENSG00000119699	transforming growth factor beta 3 Source HGNC Symbol Acc HGNC 11769
G04	UPFH1132859	ENST00000645071.1	TLR4	ENSG00000136869	toll like receptor 4 Source HGNC Symbol Acc HGNC 11850
G05	UPFH1172914	ENST00000508254.5	TLR6	ENSG00000174130	toll like receptor 6 Source HGNC Symbol Acc HGNC 16711
G06	UPFH1132978	ENST00000449264.3	TNF	ENSG00000232810	tumor necrosis factor Source HGNC Symbol Acc HGNC 11892
G07	UPFH0288248	ENST00000479933.2	TNFRSF8	ENSG00000120949	TNF receptor superfamily member 8 Source HGNC Symbol Acc HGNC 11923
G08	UPFH0607162	ENST00000615230.4	TNFRSF9	ENSG00000049249	TNF receptor superfamily member 9 Source HGNC Symbol Acc HGNC 11924
G09	UPFH0265690	ENST00000281834.4	TNFSF4	ENSG00000117586	TNF superfamily member 4 Source HGNC Symbol Acc HGNC 11934
G10	UPFH0496191	ENST00000525621.5	TYK2	ENSG00000105397	tyrosine kinase 2 Source HGNC Symbol Acc HGNC 12440
G11	UPFH0281656	ENST00000425836.6	VEGFA	ENSG00000112715	vascular endothelial growth factor A Source HGNC Symbol Acc HGNC 12680
G12	UPFH0134182	ENST00000553625.5	YY1	ENSG00000100811	YY1 transcription factor Source HGNC Symbol Acc HGNC 12856
H01	UPFH1132936	ENST00000646664.1	ACTB	ENSG00000075624	actin beta Source HGNC Symbol Acc HGNC 132
H02	UPFH1132937	ENST00000544417.5	B2M	ENSG00000166710	beta-2-microglobulin Source HGNC Symbol Acc HGNC 914
H03	UPFH1132938	ENST00000229239.10	GAPDH	ENSG00000111640	glyceraldehyde-3-phosphate dehydrogenase Source HGNC Symbol Acc HGNC 4141
H04	UPFH1132939	ENST00000298556.8	HPRT1	ENSG00000165704	hypoxanthine phosphoribosyltransferase 1 Source HGNC Symbol Acc HGNC 5157
H05	UPFH1132941	ENST00000392514.9	RPLP0	ENSG00000089157	ribosomal protein lateral stalk subunit P0 Source HGNC Symbol Acc HGNC 10371
H06	UPFH1126608	UPL_HGDC	HGDC	UPL_HGDC	Human Genomic DNA Contamination
H07	UPFH1126606	UPL_QIC	QIC	UPL_QIC	QuantiNova Internal Control
H08	UPFH1126606	UPL_QIC	QIC	UPL_QIC	QuantiNova Internal Control
H09	UPFH1126606	UPL_QIC	QIC	UPL_QIC	QuantiNova Internal Control
H10	UPFH1126605	UPL_PPC	PPC	UPL_PPC	Positive PCR Control
H11	UPFH1126605	UPL_PPC	PPC	UPL_PPC	Positive PCR Control
H12	UPFH1126605	UPL_PPC	PPC	UPL_PPC	Positive PCR Control



Related products

Product	Contents	Cat. no.
QuantiNova LNA Probe PCR QC Panel	These panels are designed to assess the quality of RNA samples before characterization using QuantiNova LNA Probe PCR Focus Panels; available in 96-well, 384-well, and Rotor-Disc 100 formats	249945
QuantiNova Reverse Transcription Kit (10)*	For 10 x 20 μ l reactions: 20 μ l 8x gDNA Removal Mix, 10 μ l Reverse Transcription Enzyme, 40 μ l Reverse Transcription Mix (containing RT primers), 20 μ l Internal Control RNA, 1.9 ml RNase-Free Water	205410
QuantiNova Probe RT-PCR Kit (100)*	For 100 x 20 μ l reactions: 1 ml QuantiNova Probe RT-PCR Master Mix, 20 μ l QuantiNova Probe RT Mix, 20 μ l Internal Control RNA, 500 μ l Yellow Template Dilution Buffer, 250 μ l ROX Reference Dye, 1.9 μ l RNase-Free Water	208352
QuantiNova Probe PCR Kit (100)*	For 100 x 20 μ l reactions: 1 ml 2x QuantiNova Probe PCR Master Mix, 500 μ l QuantiNova Yellow Template Dilution Buffer, 250 μ l QN ROX Reference Dye, 1.9 ml Water	208252

*Larger kit sizes available.

The QuantiNova LNA Probe PCR Focus Panels are intended for molecular biology applications. These products are not intended for the diagnosis, prevention or treatment of a disease.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN®, LNA®, QuantiNova®, Sample to Insight® (QIAGEN Group); SYBR® (Life Technologies Corp.). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

09/2019 © 2019 QIAGEN, all rights reserved.